



Frank Bacon Machinery Sales Company, Inc.

4433 East Eight Mile Road Warren, MI 48090-0886 (586) 756-4280 Fax (586) 756-8009 Toll Free (800) 394-4280

e-mail: sales@frankbacon.com On The Web: www.frankbacon.com

Stock #	DESCRIPTION
12240-03	<p>Buehler Phoenix 4000</p> <p>(1) USED Buehler Phoenix 4000 Sample Preparation System, Semi-Automatic or Manual variable speed Grinder/Polisher with Integral Power Head. S/N PHX-10331, Year 1995</p> <p>Buehler Phoenix 4000</p> <p>SPECIFICATIONS:</p> <p>Wheel Diameter: 12"</p> <p>Wheel Speed (RPM): 150</p> <p>Specimen Holder Speed (RPM): 50 – 600 (Both directions)</p> <p>Approximate Dimensions (W x D x H): 24" x 26" x 23"</p> <p>Approximate Weight (lbs): 210</p> <p>Electrical Requirements (V, Hz, PH): 230, 50-60, 3</p> <p>Motor (HP): 65</p> <p>Utility Requirements:</p> <p>Water: 15 – 85 PSI</p> <p>Air: 85 PSI</p> <p>Drain</p> <p>EQUIPPED WITH:</p> <ul style="list-style-type: none">• Single and Central Force Sample Operation• Touch Pad Controls with LED Display of Force, RPM and Time• Four 12" Diameter Platens• Splash Ring• Cat. # 60-5282 Central Force Specimen Holder 1.25" diameter specimens• Cat. # 60-5283 Central Force Specimen Holder 1.5" diameter specimens• User's Manual <p>FOB Warren, MI</p> <p>Cleaned, Serviced & Guaranteed</p>

Guarantee:

Unless otherwise specified, every machine is offered with the standard **MDNA** (*MACHINERY DEALERS NATIONAL ASSOCIATION*) Return privilege to ensure your complete satisfaction. If the machine is un-satisfactory it may be returned to our warehouse, freight prepaid and in the original condition within 30 Days of shipment for a full refund less the cost of SPECIALIZED EQUIPMENT, NEW ITEMS AND RE-CERTIFICATION COSTS, When applicable.

Care is taken to provide accurate specifications. However, Critical areas should be verified by Inspection.

MDNA | AMDA | MMA | ASTM | SAE | ASM